Validation of the Semiautomatic Quantification of ¹⁸F-Fluoride PET/CT Whole-Body Skeletal Tumor Burden

Ana E. Brito^{1–3}, Felipe Mourato¹, Allan Santos^{3,4}, Camila Mosci³, Celso Ramos^{3,4}, and Elba Etchebehere^{3,4}

¹Real Nuclear of Real Hospital Português de Beneficência em Pernambuco, Recife, Brazil; ²Keizo Asami Immunopathology Laboratory, Federal University of Pernambuco, Recife, Brazil; ³Division of Nuclear Medicine, Department of Radiology, University of Campinas, Campinas, Brazil; and ⁴MND Campinas, Campinas, Brazil

Our purpose was to validate a semiautomatic quantification of the skeletal tumor burden on ¹⁸F-fluoride PET/CT using manual quantification as a reference. Methods: We quantified 51 18Ffluoride PET/CT examinations performed on female breast cancer patients. Clinical information (age: time of disease presentation: presence of visceral metastases; and time to death, progression. or a bone event) was recorded. The total volume of ¹⁸F-fluorideavid skeletal metastases and the total activity of ¹⁸F-fluoride-avid metastases were calculated manually and semiautomatically. Results: Manual and semiautomatic metrics correlated strongly (P < 0.0001; 95% confidence interval, 0.9300-0.9769). On multivariable analysis, the semiautomatic measures of total activity for ¹⁸F-fluoride-avid metastasis correlated significantly with overall survival (P = 0.0001) and progression-free survival (P = 0.0001) 0.0006). Approximate times for calculating skeletal tumor burden (semiautomatic vs. manual) were, respectively, 30 s versus 321 s in patients with fewer than 5 metastases, 120 s versus 640 s in patients with 5–10 metastases, and 240 s versus 1207s in patients with more than 10 metastases. Conclusion: Semiautomatic quantification of whole-body ¹⁸F-fluoride PET/CT skeletal tumor burden can replace manual quantification in breast cancer patients and is a strong independent biomarker of prognosis.

Key Words: NaF PET/CT; breast cancer; skeletal tumor burden; semiautomatic quantification; manual quantification

J Nucl Med Technol 2018; 46:378-383

DOI: 10.2967/jnmt.118.211474

Compared with conventional bone scintigraphy or SPECT/CT, ¹⁸F-fluoride PET/CT is more sensitive and specific in the detection of osteoblastic metastases because it detects a greater number of potentially metastatic lesions (*1–3*). Because of the clear benefits of ¹⁸F-fluoride PET/CT imaging, it seems logical to take maximum advantage of this novel tracer when evaluating patients with breast and prostate cancer. For example, the determination of ¹⁸F-fluoride PET/CT whole-body skeletal tumor burden is able to objectively assess response to treatment

Received Mar. 13, 2018; revision accepted Jun. 29, 2018. For correspondence or reprints contact: Elba Etchebehere, Division of Nuclear Medicine, Department of Radiology, Rua Vital Brasil 251, Universidade

Estadual de Campinas, Campinas, SP, Brazil 13083-888. E-mail: elba.etchebehere@gmail.com

Published online Aug. 3. 2018.

COPYRIGHT © 2018 by the Society of Nuclear Medicine and Molecular Imaging.

and is a better prognostic indicator than visual analysis (4,5). Investigators have demonstrated a prognostic role for ¹⁸F-fluoride PET/CT in breast (6) and prostate (7) cancers. Therefore, quantitative parameters with ¹⁸F-fluoride PET/CT can be important for clinical evaluation of patient outcome and for clinical trials, by eliminating subjectivity and assessing response to therapy. Despite these important benefits, manual calculation of ¹⁸F-fluoride PET/CT skeletal tumor burden can be time-consuming. To routinely implement semiautomatic quantification of ¹⁸F-fluoride PET/CT skeletal tumor burden, this tool must be validated by comparison to manual quantification.

The aim of this investigation was to validate semiautomatic quantification of skeletal tumor burden on ¹⁸F-fluoride PET/CT against manual quantification.

MATERIALS AND METHODS

Patients

We reviewed all ¹⁸F-fluoride PET/CT examinations performed on 51 female patients with breast cancer and suspected bone metastases.

The local Institutional Review Board approved this retrospective study (approval 46/2016), and the requirement to obtain informed consent was waived. Clinical information (age; time of disease presentation; presence of visceral metastases; and time to death, progression, or a bone event) for each patient was also recorded.

The median follow-up time starting from the ¹⁸F-fluoride PET/CT study was 15.4 mo (range, 1.6–83 mo). The mean age of the 51 patients was 56.7 y (range, 30–82 y).

¹⁸F-Fluoride PET/CT Acquisition

The studies were performed on 2 PET/CT scanners (Biograph True-Point 64 and Biograph 16; Siemens Healthcare) 45 min after a 3.7 MBq/kg intravenous injection of ¹⁸F-sodium fluoride. PET whole-body images (from the skull to the feet) were acquired in 3-dimensional mode at a rate of 90 s/bed position and then underwent iterative reconstruction (2 iterations of 8 subsets and a gaussian filter) or TrueX (Siemens) reconstruction plus time-of-flight reconstruction (2 iterations of 21 subsets with a gaussian filter). The CT acquisition parameters included a slice thickness of 5 mm, 120 kV, or CARE Dose4D (Siemens), with no intravenous CT contrast administration.

Image Analysis and Quantification Parameters

Two experienced observers visually analyzed all images.

Two forms of volumetric quantification were undertaken in all patients to determine whole-body skeletal tumor burden: a manual quantification using syngo.via (Siemens Medical Solutions) and a semiautomatic quantification using Multi-Foci Segmentation in syngo.via.

To standardize quantification in both methods, we established an SUV_{max} cutoff of 10. We were able to calculate the skeletal tumor burden with this cutoff because it has been shown to exclude normal bone in 98% of patients (8).

Using this SUV_{max} cutoff, we calculated 2 whole-body skeletal tumor burden metrics for the manual and semiautomatic quantifications: FTV_{10} , or the total volume of ^{18}F -fluoride—avid skeletal metastases (mL), and TLF_{10} , or the total activity of ^{18}F -fluoride—avid metastases determined by the equation $FTV_{10} \times$ the average of all the SUV_{max} (g) (8).

¹⁸F-Fluoride PET/CT Manual Quantification

To perform manual quantification, we first delineated FTV_{10} as a 3-dimensional ellipsoid volume of interest (VOI) that included all metastatic ¹⁸F-fluoride–avid skeletal lesions having an SUV_{max} above 10, with an isocontour threshold set at 41% of the SUV_{max} as recommended by guidelines (9). Subsequently, this FTV_{10} was multiplied by the average SUV_{max} (g/mL) (10) to obtain TLF_{10} . The lesion with the highest SUV_{max} was also determined. Figure 1 shows an example of the manual quantification method.

The time spent on manual quantification was defined as the interval from the manual designation of the VOI surrounding the bone metastasis (including the time to annotate all necessary data) through calculation of the manual FTV and TLF ($mFTV_{10}$ and $mTLF_{10}$).

¹⁸F-Fluoride PET/CT Semiautomatic Quantification

Semiautomatic quantification with the Multi-Foci Segmentation tool in syngo.via is feasible using any SUV cutoff. After we had determined that a cutoff SUV_{max} of 10 would be used, the software automatically generated VOIs surrounding all bone lesions with an SUV greater than 10. It was important, at that point, to verify that all delineated VOIs actually corresponded to bone metastases; if any did not, we manually excluded them. Subsequently, the skeletal tumor burden was determined according to the selected tumor-related areas (Fig. 2), and the software automatically provided the TLF_{10} and FTV_{10} parameters.

The time spent on semiautomatic quantification was defined as the interval from the automatic designation of the VOI surrounding the bone metastasis (including the time to manually exclude all nontumor areas) through calculation of the semiautomatic FTV and TLF ($saFTV_{10}$ and $saTLF_{10}$).

Statistical Analysis

The Pearson test, the Bland–Altman test (11,12), and mountain plots (13) were used to test the similarity between the 2 quantification methods. The Bland–Altman test was the default method to compare the 2 measurements; mountain plots were used as a complementary tool because they better visualize the differences.

Cox proportional hazards regression was used to evaluate the relationship between the whole-body skeletal tumor burden metrics and the following variables: outcome, age, time of disease onset, highest SUV_{max} , and presence of visceral metastases.

We studied the following outcome measures: overall survival, time to progression, and time to a bone event. Overall survival was defined as the interval from the ¹⁸F-fluoride PET/CT acquisition until death. Time to progression was defined as the interval from the ¹⁸F-fluoride PET/CT acquisition until any evidence of disease progression. Time to a bone event was defined as the interval from the ¹⁸F-fluoride PET/CT acquisition until the occurrence of any




FIGURE 1. Example of manual quantification for one patient, with $mFTV_{10}$ corresponding to sum of all VOIs (81.38 mL) and $mTLF_{10}$ corresponding to $mFTV_{10}$ multiplied by average of all SUV_{max} (1,266.27 g).




FIGURE 2. (A) After threshold selection (SUV_{max} of 10), several VOIs are automatically created. (B) On visual inspection, VOIs not related to metastases are deleted. In this patient, as indicated by white arrow in coronal section on CT scan, uptake was noted in calcified portion of liver, which was a sequela of radioablation (L11VOI1); this lesion was manually excluded. (C) Patient's final $saTLF_{10}$ was 3,969 g.

bone event. Even though bone events indicate progression of the disease, these events were calculated separately because of morbidity and mortality.

All statistical analyses were performed using MedCalc software. The level of significance was set at 5%.

RESULTS

During follow-up, 21 patients died (41.2%) and 29 progressed (56.9%). Nine patients who progressed had bone events (17.6%). One patient had no progression and remained alive throughout the follow-up.

¹⁸F-Fluoride PET/CT Manual and Semiautomatic Quantifications

Because of the strong correlation between $mTLF_{10}$ and $mFTV_{10}$ ($\rho=0.8117$; P<0.0001; 95% confidence interval, 0.6905–0.8885), only $mTLF_{10}$ was used for subsequent analyses (Fig. 3). The same was true for $saTLF_{10}$ and $saFTV_{10}$ ($\rho=0.9234$; P<0.0001; 95% confidence interval, 0.8690–0.9558), with only $saTLF_{10}$ being used for subsequent analyses (Fig. 4).

The median $mTLF_{10}$ was 255 g (range, 18–23.027 g), and the median $saTLF_{10}$ was 574 g (range, 8–32.225 g).

The Pearson correlation coefficient between $mTLF_{10}$ and $saTLF_{10}$ was 0.9596 (P < 0.0001; 95% confidence interval, 0.9300–0.9769). The Bland–Altman test showed similarity between the manual and semiautomatic methods. There was only one outlier, as seen in Figure 5. That patient had a high tumor burden on semiautomatic quantification and a low tumor burden on manual quantification (TLF_{10} , 4,523 and 327.2 g, respectively). The mountain plots showed the center to be near zero, with a left tail (Fig. 6).

Because of the strong correlation between $mTLF_{10}$ and $saTLF_{10}$, a multivariable analysis was undertaken solely on the semiautomatic measures. The analysis showed that $saTLF_{10}$ correlated significantly with overall and progression-free survival (Table 1). None of the variables were

associated with bone events. The highest SUV_{max} did not bear a significant prognostic value.

Approximate times for quantifying skeletal tumor burden (semiautomatic vs. manual) were, respectively, 30 s versus 321 s in patients with fewer than 5 metastases, 120 s versus 640 s in patients with 5–10 metastases, and 240 s versus 1,200 s in patients with more than 10 lesions.

DISCUSSION

Semiautomatic quantification of wholebody skeletal tumor burden with ¹⁸F-fluoride PET/CT (TLF₁₀ and FTV₁₀) has been shown to be a strong and independent prognostic biomarker in breast (6) and prostate (7,14) cancers. To the best

of our knowledge, semiautomatic quantification of ¹⁸F-fluoride PET/CT had not previously been validated for clinical use. Therefore, we manually quantified whole-body skeletal tumor burden (the method of quantification most often used in clinical practice) and compared it with semiautomatic quantification to validate the latter method.

Like studies comparing manual with semiautomatic quantification for ¹⁸F-FDG PET/CT—studies that revealed a high level of agreement (15)—our study for ¹⁸F-fluoride PET/CT obtained comparable results. In our study, mountain plots were centered over zero for both the manual and the semiautomatic methods, with a small left tail, showing that although the methods had some differences, they correlated strongly overall. These differences were secondary to an underestimation of the manually quantified values in a few cases, when compared with the semiautomatic method. Furthermore, the visual assessment of the Bland–Altman scatterplots




FIGURE 3. Graph showing strong correlation between values for mTLF₁₀ and mFTV₁₀ ($\rho = 0.8117; P < 0.0001; 95\%$ confidence interval, 0.6905–0.8885).




FIGURE 4. Graph showing strong correlation between values for $saTLF_{10}$ and $saFTV_{10}$ (p = 0.9234; P < 0.0001; 95% confidence interval, 0.8690–0.9558).

showed uniform bias, and we can thus conclude that these metrics are similar, equivalent, and strongly correlated.

Although both methods of quantification are related, the manual method is time-consuming. In addition, the manually drawn VOIs are elliptic; in patients with multiple lesions, the observer might reduce the size of these VOIs to avoid overlap of VOIs, and this reduction might reduce the total volume of lesions. That is why there was such a discrepancy between the median $mTLF_{10}$ (255 g) and the median $saTLF_{10}$ (574 g). This issue does not occur in the semiautomatic quantification because the software draws the VOIs according to the exact shape of the lesions, without overlap.

The semiautomatic quantification, therefore, clearly becomes the method of choice for whole-body tumor determination. Furthermore, it was necessary to confirm that $saTLF_{10}$ and $saFTV_{10}$ were prognostic factors in breast cancer. Specifi-




FIGURE 5. Bland–Altman test showing similarity between manual and semiautomatic quantifications.




FIGURE 6. Mountain plot showing similarity between manual and semiautomatic quantifications, with peak near zero and small left tail.

cally, when quantifying osteoblastic bone metastases on the ¹⁸F-fluoride PET/CT studies, we applied a cutoff SUV_{max} of 10, which has been shown to separate normal bone from malignant lesions in approximately 98% of patients (8,16). This cutoff is not undisputed, and some authors prefer a cutoff of 15 to minimize nonpathologic and normal-bone uptake (17,18). There are no comparative studies between the 2 cutoffs (10 and 15) to decide which is better. However, we believe that increasing the cutoff reduces the sensitivity (while increasing the specificity), potentially leading to erroneous measurement of skeletal tumor burden. After semiautomatic quantification, a visual inspection is warranted to manually exclude VOIs placed in areas unrelated to bone metastases. Hence, increasing the SUV_{max} cutoff (to 15, for instance) may render this task more time-consuming because not only will areas unrelated to bone metastases have to be excluded but areas related to bone metastases will have to be included.

One limitation of both quantification methods relates to metastatic lesions whose activity is near the level of physiologic activity. For instance, in a metastasis in the

TABLE 1
Correlation of Clinical and Imaging Variables to Overall
Survival and Time to Progression on Multivariable Analysis

Variable	95%CI	P
Overall survival		
saTLF ₁₀	1.0003-1.0009	0.0001
Visceral metastases	1.8334-24.6637	0.0041
Age	1.0042-1.1148	0.0342
Time to progression		
saTLF ₁₀	1.0001-1.0002	0.0006
Age	0.9359-0.9949	0.0342
-		

CI = confidence interval.

pubic bone, which is near the bladder, a perfect delimitation was sometimes not possible because of the physiologic urinary excretion. Manual quantification also has limitations, in addition to being laborious. First, because it is observer-dependent, the VOI may vary widely, potentially influencing the whole-body quantification results. Second, the number of quantifiable lesions is not perfect. Although manual quantification was attempted in all patients regardless of the number of lesions, in some patients with a very high tumor burden (above 20 VOIs delineated), the software processing began to slow (most likely because of work overload) and the workstation had to be restarted. This problem did not occur with the semiautomatic quantification method; therefore, that method has the additional benefit of being able to include all patients, even those with a high volume of disease. A recent study showed that manual volumetric quantification of ⁶⁸Ga-PSMA PET/CT correlated with serum prostate-specific antigen values in patients with suspected biochemical recurrence of prostate cancer (19), but unfortunately, the authors excluded from the analysis patients with more than 10 lesions, claiming that manual measurement of tumor volumes was time-consuming and labor-intensive. To date, there is only one clear advantage of manual over semiautomatic quantification: the manual method allows drawing of an irregular ROI. In that scenario, semiautomatic quantification software still needs improvement, as there is no freehand drawing tool for VOIs, allowing only spheric VOIs and their variations. Segmentation is a necessary step in semiautomatic quantification and requires human interaction. Totally automatic forms of volumetric quantification will dispense with the segmentation process; such procedures are under study and will require extensive validation before clinical use (20).

One clear strength of this study was the possibility of addressing the clinical relevance of these quantitative findings. In our retrospective analysis of female breast cancer patients, semiquantitative determination of whole-body tumor burden on ¹⁸F-fluoride PET/CT independently correlated with overall survival (6). This finding demonstrates the potential of this approach for prognostication of these patients. Although we have validated this approach for breast cancer patients, our findings need confirmation for other cancer types.

Guidelines recommend that tumor burden be reported when available (9), but such reporting is feasible only for localized solid tumors. In current clinical practice, information on $^{18}\mbox{F-fluoride PET/CT}$ skeletal tumor burden cannot be reported for patients with widespread disease such as lymphomas and bone metastases from breast and prostate cancer and must be reserved for research. A quantitative metric widely used in clinical practice because of its ease of measurement is SUV $_{\rm max}$. In our analyses, the highest SUV $_{\rm max}$ on $^{18}\mbox{F-fluoride PET/CT}$ examination was not a significant prognostic factor, as also shown by other investigators (6), even though the SUV $_{\rm max}$ of the primary lesion is a

prognostic factor in breast cancer patients on ^{18}F -FDG PET/CT (21,22). The lack of correlation between the highest SUV_{max} and outcome justifies calculation of whole-body skeletal tumor burden and not just the highest SUV_{max} .

In daily clinical practice, quantification of whole-body skeletal tumor burden has to be in at least a semi-automatic form. In this respect, the quantification proposed here is feasible, fast, and reproducible. These advantages open the door to new clinical possibilities in the theranostics setting (7,14) and in evaluating response to therapy (17).

CONCLUSION

Semiautomatic determination of whole-body skeletal tumor burden might be able to replace manual quantification for ¹⁸F-fluoride PET/CT in breast cancer patients. Furthermore, it is a strong independent prognostic imaging biomarker. These findings need confirmation in other cancer types.

DISCLOSURE

No potential conflict of interest relevant to this article was reported.

REFERENCES

- Bortot DC, Amorim BJ, Oki GC, et al. ¹⁸F-fluoride PET/CT is highly effective for excluding bone metastases even in patients with equivocal bone scintigraphy. *Eur J Nucl Med Mol Imaging*. 2012;39:1730–1736.
- Ota N, Kato K, Iwano S, et al. Comparison of ¹⁸F-fluoride PET/CT, ¹⁸F-FDG PET/CT and bone scintigraphy (planar and SPECT) in detection of bone metastases of differentiated thyroid cancer: a pilot study. *Br J Radiol*. 2014; 87:20130444.
- Even-Sapir E, Metser U, Mishani E, Lievshitz G, Lerman H, Leibovitch I.
 The detection of bone metastases in patients with high-risk prostate cancer: 99mTc-MDP planar bone scintigraphy, single- and multi-field-of-view SPECT, 18F-fluoride PET, and 18F-fluoride PET/CT. J. Nucl. Med. 2006;47:287–297.
- Juweid ME, Wiseman GA, Vose JM, et al. Response assessment of aggressive non-Hodgkin's lymphoma by integrated International Workshop Criteria and fluorine-18-fluorodeoxyglucose positron emission tomography. *J Clin Oncol*. 2005;23:4652–4661.
- Lin C, Itti E, Haioun C, et al. Early ¹⁸F-FDG PET for prediction of prognosis in patients with diffuse large B-cell lymphoma: SUV-based assessment versus visual analysis. *J Nucl Med.* 2007;48:1626–1632.
- Brito AE, Santos A, Sasse AD, et al. ¹⁸F-fluoride PET/CT tumor burden quantification predicts survival in breast cancer. *Oncotarget*. 201730;8:36001–36011.
- Etchebehere EC, Araujo JC, Fox PS, Swanston NM, Macapinlac HA, Rohren EM. Prognostic factors in patients treated with ²²³Ra: the role of skeletal tumor burden on baseline ¹⁸F-fluoride PET/CT in predicting overall survival. *J Nucl Med.* 2015;56:1177–1184.
- Rohren EM, Etchebehere EC, Araujo JC, et al. Determination of skeletal tumor burden on ¹⁸F-fluoride PET/CT. J Nucl Med. 2015;56:1507–1512.
- Boellaard R, Delgado-Bolton R, Oyen WJG, et al. FDG PET/CT: EANM procedure guidelines for tumour imaging: version 2.0. Eur J Nucl Med Mol Imaging. 2015;42:328–354.
- Kinahan PE, Fletcher JW. Positron emission tomography-computed tomography standardized uptake values in clinical practice and assessing response to therapy. Semin Ultrasound CT MR. 2010;31:496–505.
- Bland JM, Altman DG. Statistical methods for assessing agreement between two methods of clinical measurement. *Lancet*. 1986;1:307–310.
- Bland JM, Altman DG. Measuring agreement in method comparison studies. Stat Methods Med Res. 1999;8:135–160.
- Krouwer JS, Monti KL. A simple, graphical method to evaluate laboratory assays. Eur J Clin Chem Clin Biochem. 1995;33:525–527.
- Fosbøl MØ, Petersen PM, Kjaer A, Mortensen J. Radium-223 therapy of advanced metastatic castration-resistant prostate cancer: quantitative assessment of

- skeletal tumor burden for prognostication of clinical outcome and hematological toxicity. *J Nucl Med.* 2018;59:596–602.
- Torigian DA, Lopez RF, Alapati S, et al. Feasibility and performance of novel software to quantify metabolically active volumes and 3D partial volume corrected SUV and metabolic volumetric products of spinal bone marrow metastases on ¹⁸F-FDG-PET/CT. Hell J Nucl Med. 2011;14:8–14.
- Lapa P, Marques M, Costa G, Iagaru A, Pedroso de Lima J. Assessment of skeletal tumour burden on ¹⁸F-NaF PET/CT using a new quantitative method. Nucl Med Commun. 2017;38:325–332.
- Harmon SA, Perk T, Lin C, et al. Quantitative assessment of early [18F]sodium fluoride positron emission tomography/computed tomography response to treatment in men with metastatic prostate cancer to bone. *J Clin Oncol*. 2017;35:2829–2837.
- Lin C, Bradshaw T, Perk T, et al. Repeatability of quantitative ¹⁸F-NaF PET: a multicenter study. J Nucl Med. 2016;57:1872–1879.

- Schmuck S, von Klot CA, Henkenberens C, et al. Initial experience with volumetric ⁶⁸Ga-PSMA I&T PET/CT for Assessment of whole-body tumor burden as a quantitative imaging biomarker in patients with prostate cancer. *J Nucl Med*. 2017;58:1962–1968.
- Taghanaki SA, Duggan N, Ma H, et al. Segmentation-free direct tumor volume and metabolic activity estimation from PET scans. Comput Med Imaging Graph. 2018:63:52–66
- Yoo J, Yoon H-J, Kim BS. Prognostic value of primary tumor SUVmax on F-18 FDG PET/CT compared with semi-quantitative tumor uptake on Tc-99m sestamibi breastspecific gamma imaging in invasive ductal breast cancer. Ann Nucl Med. 2017;31:19–28.
- García García-Esquinas M, García-Sáenz JA, Arrazola García J, et al. ¹⁸F-FDG PET-CT imaging in the neoadjuvant setting for stages II-III breast cancer: association of loco-regional SUVmax with classical prognostic factors. *Q J Nucl Med Mol Imaging*. 2014;58:66–73.

Erratum

The article "Exposure to Technologists from Preparing and Administering Therapeutic ¹³¹I: How Frequently Should We Bioassay?" by Kopisch et al. (*J Nucl Med Technol.* 2011;39:60–62) contained a mathematical error that was addressed in a June 2011 Erratum. However, additional corrections in the article are required based on the mathematical error correction in the "Results" section. The corrected paragraphs appear below:

Results

From these data, one learns that an average air concentration of **2.4E–06 kBq/mL** (**6.4E–08 \muCi/mL**) of air can be expected from the handling and administration of a dose of 5.74 GBq (155 mCi) of ¹³¹I. The NRC assumption for its derived air concentration calculations is that an average worker inhales approximately 20 L of air per minute. A technologist utilizing a full 10 min of ¹³¹I handling for a procedure would inhale about **200 L of air**. One could project a total ingestion for the technologist of **0.481 kBq/mL** (**0.013 \muCi**) during such a procedure.

Discussion

First paragraph: Table 2 summarizes the number of procedures and the number of participating technologists at each of the study locations. The average number of 131 I procedures performed by each technologist in this study was 4. Therefore, the average 131 I dose a technologist in this study received in a year was 4×0.481 kBq (0.013 μ Ci), or 19.2 kBq (0.052 μ Ci), well below the NRC monitoring guideline of 185 kBq (5 μ Ci) per year. The actual dose would probably be lower because this estimate assumes an average dose activity of 5.74 GBq (155 mCi).

Third paragraph: This small-scale study's results indicate a typical annual intake that is well below the level that the NRC advises as a trigger level for bioassay monitoring. The study results indicate that one would have to administer close to **2220 GBq (60,000 mCi)** of ¹³¹I in 1 y to reach the NRC trigger limit for bioassay.

Conclusion

First paragraph: This small study showed an average 131 I inhalation intake of **0.481 kBq (0.013 \muCi)** during administration of 5.74 GBq (155 mCi) of 131 I in capsule form. This value allows for a full 10 min to handle, assay, and administer the dose. On the basis of this small-scale study, a technologist would have to administer approximately **185 GBq (5,000 mCi)** of 131 I every month to trigger the NRC threshold of 10% ALI for 131 I. This is the trigger level the NRC recommends for bioassay of occupational workers.

In addition, the last column in **Table 1** requires corrections and should read:

¹³¹ I air concentration (Bq/mL)
6.73E-05
8.18E-05
6.03E-05
2.90E-03
1.05E-02
8.23E-04
2.13E-03